

## Product Datasheet

### Bet v 1.0101 protein orb424316

<b>Description</b>	Recombinant of Bet v 1.0101 protein
<b>Conjugation</b>	Unconjugated
<b>Target</b>	Bet v 1.0101
<b>Form/Appearance</b>	Sterile Filtered clear solution. Bet v 1.0101 protein solution is supplied in 20mM HEPES buffer pH-8.0, 200mM NaCl and 20% glycerol.
<b>Storage</b>	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
<b>Note</b>	For research use only.
<b>Purity</b>	>95.0% as determined by SDS-PAGE.
<b>Source</b>	Sf9 insect cells.
<b>Expiration Date</b>	6 months from date of receipt.