

Product Datasheet

Borrelia DbpB protein orb424373

Description	Recombinant of Borrelia DbpB protein
Conjugation	Unconjugated
Tested Applications	WB
Target	Borrelia DbpB
Form/Appearance	Sterile Filtered clear solution. Borrelia DbpB (1.11mg/1ml) is supplied in 16mM HEPES buffer pH-8.0, 300mM NaCl and 20% glycerol.
Storage	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
Note	For research use only.
Purity	>80.0% as determined by SDS-PAGE.
Source	Escherichia Coli.
Expiration Date	6 months from date of receipt.