

Product Datasheet

Human CAPS protein orb424501

Description	Recombinant of human CAPS protein
Reactivity	Human
Conjugation	Unconjugated
Target	CAPS
Form/Appearance	Sterile Filtered colorless solution. The CAPS solution (1 mg/ml) contains 20mM Tris-HCl buffer (pH8.0), 10% glycerol, 2mM DTT and 100mM NaCl.
Storage	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
Note	For research use only.
Protein Sequence	MGSSHHHHHH SSGLVPRGSH MDAVDATMEK LRAQCLSRGA SGIQGLARFF RQLDRDGSRS LDADEFROGL AKLGLVLDQA EAEGVCRKWD RNGSGTLDLE EFLRALRPPM SQAREAVIAA AFAKLDRSGD GVVTVDLLRG VYSGRAHPKV RSGEWTEDEV LRRFLDNFDS SEKDGQVTLA EFQDYYSQVVS ASMNTDEEFV AMMTSAWQL.
Purity	>95.0% as determined by SDS-PAGE.
Source	Escherichia Coli.
Expiration Date	6 months from date of receipt.