

## Product Datasheet

### Human CD46 (Sf9) protein orb424611

<b>Description</b>	Recombinant of human CD46 (Sf9) protein
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Target</b>	CD46 (Sf9)
<b>Form/Appearance</b>	Sterile filtered colorless solution. CD46 protein solution (0.5mg/ml) contains Phosphate buffered saline (pH7.4) and 10% glycerol.
<b>Storage</b>	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
<b>Note</b>	For research use only.
<b>Protein Sequence</b>	ADLCEEPPTF EAMELIGKPK PYEIGERVD YKCKKGYFYI PPLATHTICD RNHTWLPVSD DACYRETCPY IRDPLNGQAV PANGTYEFGY QMHFICNEGY YLIGEEILYC ELKGSVAIWS GKPPICEKVL CTPPKIKNG KHTFSEVEVF EYLDAVTYSC DPAPGDPFVS LIGESTIYCG DNSVWSRAAP ECKVVKCRFP VVENGKQISG FGKKFYKAT VMFECDKGFY LDGSDTIVCD SNSTWDPPVP KCLKGPRPTY KPPVSNYPGY PKPEEGILDS LDHHHHHH
<b>Purity</b>	>90.0% as determined by SDS-PAGE.
<b>Source</b>	Sf9, Baculovirus cells.
<b>Expiration Date</b>	6 months from date of receipt.