

## Product Datasheet

### Human CDK3 protein orb424671

<b>Description</b>	Recombinant of human CDK3 protein
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Target</b>	CDK3
<b>Form/Appearance</b>	Sterile Filtered colorless solution. The CDK3 solution (0.25mg/ml) contains 20mM Tris-HCl buffer (pH 8.0), 1M Urea and 10% glycerol.
<b>Storage</b>	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
<b>Note</b>	For research use only.
<b>Protein Sequence</b>	MGSSHHHHHH SSGLVPRGSH MGSMDMFQKV EKIGEGTYGV VYKAKNRETG QLVALKKIRL DLEMEGVPST AIREISLLKE LKHPNIVRLL DVVHNERKLY LVFEFLSQDL KKYMDSTPGS ELPLHLIKSY LFQLLQGVSF CHSHRVIHRD LKPQNLLINE LGAIKLADFG LARAFGVPLR TYTHEVVTW YRAPEILLGS KFYTTAVDIW SIGCIFAEMV TRKALFPGDS EIDQLFRIFR MLGTPSEDW PGVTQLPDYK GSFPKWTRKG LEEIVPNLEP EGRDLLMQLL QYDPSQRITA KTALAHYPYFS SPEPSAARQ YVLQRFRRH.
<b>Purity</b>	>85% as determined by SDS-PAGE.
<b>Source</b>	Escherichia Coli.
<b>Expiration Date</b>	6 months from date of receipt.