

Product Datasheet

Human CEACAM7 protein orb424692

Description	Recombinant of human CEACAM7 protein
Reactivity	Human
Conjugation	Unconjugated
Target	CEACAM7
Form/Appearance	Sterile Filtered clear solution. CEACAM7 protein solution (0.25mg/ml) contains 20mM Tris-HCl buffer (pH 8.0), 0.15M NaCl and 10% glycerol.
Storage	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
Note	For research use only.
Protein Sequence	MGSSHHHHHH SSGLVPRGSH MGSTNIDVVP FNVAEGKEVL LVVHNESQNL YGYNWYKGER VHANYRIIGY VKNISQENAP GPAHNGRETI YPNGTLLIQN VTHNDAGIYT LHVIKENLVN EEVTRQFYVF SEPPKPSITS NNFNPVENKD IVVLTCPET QNTTYLWWVN NQSLVSPRL LLSTDNRTL VLLSATKNDIG PYECEIQNPV GASRSDPVTL NVRYESVQAS.
Purity	>85.0% as determined by SDS-PAGE.
Source	Escherichia Coli.
Expiration Date	6 months from date of receipt.