

Product Datasheet

Baby Rabbit Complement orb436707

Description	Baby Rabbit Complement
Reactivity	Rabbit
Tested Applications	Functional Assay
Hazard Information	Non-Toxic
Target	Baby Rabbit Complement
Form/Appearance	Lyophilised
Storage	Store at 4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Note	For research use only.
Purity	Baby rabbit serum
Expiration Date	12 months from date of receipt.