

Product Datasheet

E. coli IDNK protein orb426348

Description	Recombinant of e. coli IDNK protein
Reactivity	E. coli
Conjugation	Unconjugated
Target	IDNK
Form/Appearance	Sterile Filtered colorless solution. IDNK protein solution (1mg/ml) containing 20mM Tris-HCl buffer (pH 8.0), 0.15M NaCl, 10% glycerol and 1mM DTT.
Storage	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
Note	For research use only.
Protein Sequence	MGSSHHHHHH SSGLVPRGSH MGSMAGESFI LMGVSGSGKT LIGSKVAALL SAKFIDGDDL HPAKNIDKMS QGIPLSDEDR LPWLERLNDA SYSLYKKNET GFIVCSSLKK QYRDILRKGS PHVFHLWLDG DYETILARMQ RRAGHFMPVA LLKSQFEALE RPQADEQDIV RIDINHDIAN VTEQCRQAVL AIRQNRICAK EGSASDQRCE.
Purity	>95% as determined by SDS-PAGE.
Source	E.coli.
Expiration Date	6 months from date of receipt.