

Product Datasheet

Human IFNA14 protein orb426374

Description	Recombinant of human IFNA14 protein
Reactivity	Human
Conjugation	Unconjugated
Target	IFNA14
Form/Appearance	Sterile filtered colorless solution. IFNA14 protein solution (0.5mg/ml) containing 20mM Tris-HCl buffer (pH8.0), 50% glycerol, 0.2M NaCl and 2mM DTT.
Storage	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
Note	For research use only.
Protein Sequence	MGSSHHHHHH SSGLVPRGSH MGSHEMCNLSQ THSLNNRRTL MLMAQMRRIS PFCLKDRHD FEFQEEFDG NQFQKAQAIS VLHEMMQQTFLNFSTKNSSA AWDETLLLEKF YIELFQQMND LEACVIEVGVVEETPLMNE SILAVKKYFQ RITLYLMEKK YSPCAWEVVR AEIMRSLSFS TNLQKRLRRK D.
Purity	>90% as determined by SDS-PAGE.
Source	Escherichia Coli.
Expiration Date	6 months from date of receipt.