

Product Datasheet

Human IFT20 protein orb426383

Description	Recombinant of human IFT20 protein
Reactivity	Human
Conjugation	Unconjugated
Target	IFT20
Form/Appearance	Sterile Filtered colorless solution. The IFT20 solution contains 20mM Tris-HCl buffer (pH 8.0), 0.15M NaCl, 1mM DTT and 30% glycerol.
Storage	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
Note	For research use only.
Protein Sequence	MGSSHHHHHH SSGLVPRGSH MAKDILGEAG LHFDELNKLRL VLDPEVTQQT IELKEECKDF VDKIGQFQKI VGGLIELVDQ LAKEAENEKM KAIGARNLLK SIAKQREAQQ QQLQALIAEK KMQLERYRVE YEALCKVEAE QNEFIDQFIF QK
Purity	>95% as determined by SDS-PAGE.
Source	E.coli.
Expiration Date	6 months from date of receipt.