

Product Datasheet

Human IL 12 p35 protein orb426433

| | |
|------------------------|--|
| Description | Recombinant of human IL 12 p35 protein |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Target | IL 12 p35 |
| Form/Appearance | Sterile Filtered clear solution. Interleukin-12 p35 His-Tag is supplied in 1xPBS and 50% glycerol. |
| Storage | Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles. |
| Note | For research use only. |
| Purity | >95.0% as determined by SDS-PAGE. |
| Source | Escherichia Coli. |
| Expiration Date | 6 months from date of receipt. |