

Product Datasheet

Human IL 2 (His) protein orb426481

Description	Recombinant of human IL 2 (His) protein
Reactivity	Human
Conjugation	Unconjugated
Target	IL 2 (His)
Form/Appearance	Sterile Filtered clear solution. Interleukin-2 His-Tag protein is supplied in 20mM Tris-HCl pH-8 and 50% glycerol.
Storage	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
Note	For research use only.
Purity	>90.0% as determined by SDS-PAGE.
Source	Escherichia Coli.
Expiration Date	6 months from date of receipt.