

## Product Datasheet

### IL11RA protein orb426551

<b>Description</b>	Recombinant of IL11RA protein
<b>Conjugation</b>	Unconjugated
<b>Target</b>	IL11RA
<b>Form/Appearance</b>	Sterile filtered colorless solution. IL11RA protein solution (0.25mg/ml) contains Phosphate Buffered Saline (pH7.4) and 10% glycerol.
<b>Storage</b>	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
<b>Note</b>	For research use only.
<b>Protein Sequence</b>	TPCPQAWGPP GVQYGQPGRP VMLCCPGVNA GTPVSWFRDG DSRLLQGPDS GLGHRVLVAQ VDSRDEGTYV CRTLDGVFVG MVTLLKLGSP ARPEVSCQAV DYENFSCTWS PGRVSGLPTR YLTSYRKKTL PGAESQRESP STGPWPCPD PLEASRCVH GAEFWSEYRI NVTEVNPLGA STCLLDVRLQ RILRPDPPQ LRVESVPGYP RRLHASWTYP ASWRRQPHFL LKFRLQYRPA QHPAWSTVEP IGLEELITDA VAGLPHAVRV SARDFLDAGT WSAWSPEAWG TPSTGPLRDE VPDGSRGHEQ KLEAAAQEDS PAPPSPSLQP DPRPLDHRDP LEQVAVLAVE HHHHHH.
<b>Purity</b>	>85.0% as determined by SDS-PAGE.
<b>Source</b>	Sf9, Baculovirus cells.
<b>Expiration Date</b>	6 months from date of receipt.