

Product Datasheet

Mouse IL34 protein orb426595

Description	Recombinant of mouse IL34 protein
Reactivity	Mouse
Conjugation	Unconjugated
Target	IL34
Form/Appearance	Sterile Filtered colorless solution. IL34 protein solution (0.5mg/ml) in 20mM Tris-HCl buffer (pH 8.0) and 10% Glycerol.
Storage	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
Note	For research use only.
Protein Sequence	MGSSHHHHHH SSGLVPRGSH MNENLEIWTL TQDKECDLTG YLRGKLQYKN RLQYMKHYFP INYRIAVPYE GVLRVANITR LQKAHVSERE LRYLWVLVSL NATESVMDVL LEGHPSWKYL QEVQTLLENV QRSLMDVEIG PHVEAVLSLL STPGLSLKLV RPKALLDNCF RVMELLYCSC CKQSPILKWQ DCELPRLHPH SPGSLMQCTA TNVYPLSRQT PTLPGSPSS SHGSLP.
Purity	>85% as determined by SDS-PAGE.
Source	E.coli.
Expiration Date	6 months from date of receipt.