

Product Datasheet

Mouse IL36B protein orb426604

Description	Recombinant of mouse IL36B protein
Reactivity	Mouse
Conjugation	Unconjugated
Target	IL36B
Form/Appearance	Sterile Filtered White Lyophilized (freeze-dried) powder. Lyophilized from a 0.2µm filtered concentrated solution in 1
Storage	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
Note	For research use only.
Protein Sequence	MMAFPQSCV HVLPPKSIQM WEPNHNTMHG SSQSPRNYRV HDSQQMVWVL TGNTLTAVPA SNNVKPVILS LIACRDTEFQ DVKKGNLVFL GIKNRNLCFC CVEMEGKPTL QLKEVDIMNL YKERKAQKAF LFYHGIEGST SVFQSVLYPG WFIATSSIER QTIILTHQRG KLVNTNFYIE SEK
Purity	>95.0% as determined by SDS-PAGE.
Source	Escherichia Coli.
Activity	Measured by its binding ability in a functional ELISA to bind recombinant mouse IL-1 Rrp2.
Solubility (25°C)	It is recommended to reconstitute the lyophilized IL36B in sterile 18M-cm H2O not less than 100µg/ml, which can then be further diluted to other aqueous solutions.
Expiration Date	6 months from date of receipt.