

Product Datasheet

Human IL4R protein orb426614

Description	Recombinant of human IL4R protein
Reactivity	Human
Conjugation	Unconjugated
Target	IL4R
Form/Appearance	Sterile filtered colorless solution. IL4R protein solution (0.25mg/ml) contains Phosphate buffered saline (pH7.4).
Storage	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
Note	For research use only.
Protein Sequence	MKVLQEPTCV SDYMSISTCE WKMNGPTNCS TELRLLYQLV FLLSEAHTCI PENNGGAGCV CHLLMDDVVS ADNYTLDLWA GQQLLWKGSF KPSEHVKPRG PGNLTVHTNV SDTLLLTWSN PYPPDNYLYN HLTYAVNIWS ENDPADFRIY NVTYLEPSLR IAASTLKSGI SYRARVRAWAQCYNNTTWESEW SPSTKWHNSY REPFEQHLEH HHHHH
Purity	>95.0% as determined by SDS-PAGE.
Source	Sf9, Baculovirus cells.
Expiration Date	6 months from date of receipt.