

Product Datasheet

Mouse IL-6 (His) protein orb426618

Description	Recombinant of mouse IL-6 (His) protein
Reactivity	Mouse
Conjugation	Unconjugated
Target	IL-6 (His)
Form/Appearance	Sterile Filtered colorless solution. IL6 Mouse protein solution (0.25mg/ml) containing 20mM Tris-HCl (pH 8.0) and 10% glycerol.
Storage	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
Note	For research use only.
Protein Sequence	MGSSHHHHHH SSGLVPRGSH MFPTSQVRRG DFTEDTTPNR PVYTTSQVGG LITHVLWEIV EMRKELCNGN SDCMNNDDAL AENNLKLPEI QRNDGICYQTG YNQEICLLKI SSGLLEYHSY LEYMKNNLKD NKKDKARVLQ RDTETLIHIF NQEVKDLHKI VLPTPISNAL LTDKLESQKE WLRTKTIQFI LKSLEEFLLKV TLRSTRQT.
Purity	>80.0% as determined by SDS-PAGE.
Source	Escherichia Coli.
Expiration Date	6 months from date of receipt.