

## Product Datasheet

### Influenza-B Tokio protein orb426641

<b>Description</b>	Influenza-B Tokio protein
<b>Conjugation</b>	Unconjugated
<b>Target</b>	Influenza-B Tokio
<b>Form/Appearance</b>	Sterile Filtered colorless solution. The B/Tokio/53/99solution (1.7mg/ml) contains STE, 0.1% sodium azide (NaN <sub>3</sub> ) and 0.005% thimerosal.
<b>Storage</b>	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
<b>Note</b>	For research use only.
<b>Purity</b>	>90.0% as determined by Analysis by SDS-PAGE.
<b>Expiration Date</b>	6 months from date of receipt.