

Product Datasheet

Human IRF 5 protein orb426665

Description	Recombinant of human IRF 5 protein
Reactivity	Human
Conjugation	Unconjugated
Target	IRF 5
Form/Appearance	Sterile Filtered clear solution. IRF5 protein solution (0.25mg/ml) containing 20mM Tris-HCl buffer (pH 8.0), 0.15M NaCl, 30% glycerol and 1mM DTT.
Storage	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
Note	For research use only.
Protein Sequence	MRGSHHHHHH GMASMTGGQQ MGRDLYDDDD KDRWGSPTL QPPTLQPPVV LGPPAPDPSP LAPPPGNPAG FRELLSEVLE PGPLPASLPP AGEQLLPDLL I
Purity	>85% as determined by SDS-PAGE.
Source	Escherichia Coli.
Expiration Date	6 months from date of receipt.