

Product Datasheet

Human ISOC2 protein orb426669

Description	Recombinant of human ISOC2 protein
Reactivity	Human
Conjugation	Unconjugated
Target	ISOC2
Form/Appearance	Sterile Filtered colorless solution. The ISOC2 solution (0.5mg/ml) contains 20mM Tris-HCl buffer (pH 8.0), 0.4M Urea and 10% glycerol.
Storage	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
Note	For research use only.
Protein Sequence	MGSSHHHHHH SSGLVPRGSH MGSEQYPQGL GPTVPELGTE GLRPLAKTCF SMVPALQEL DSRPQLRSVL LCGIEAQACI LNTTLDLLDR GLQVHVVDVDA CSSRSQVDRL VALARMRQSG AFLSTSEGLI LQLVGDAVHP QFKEIQKLIK EPAPDSGLLG LFQGQNSLLH.
Purity	>90.0% as determined by SDS-PAGE.
Source	Escherichia Coli.
Expiration Date	6 months from date of receipt.