

Product Datasheet

Human ITGB1 (Sf9) protein orb426671

Description	Recombinant of human ITGB1 (Sf9) protein
Reactivity	Human
Conjugation	Unconjugated
Target	ITGB1 (Sf9)
Form/Appearance	Sterile Filtered colorless solution. ITGB1 protein solution (1.0mg/ml) contains Phosphate Buffered Saline (pH 7.4) and 10% glycerol.
Storage	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
Note	For research use only.
Protein Sequence	QTDENRCLKA NAKSCGECIQ AGPNCGWCTN STFLQEGMPT SARCDLLEAL KKKGCPPDDI ENPRGSKDIK KNKNVTNRSK GTAEKLPED ITQIQPQLV LRLRSGEPQT FTLKFKRAED YPIDLYLMD LSYSMKDDLE NVKSLGTDLM NEMRRITSDF RIGFGSFVEK TVMPYISTTP AKLRNPCTSE QNCTSPFSYK NVLSLTNKGE VFNELVGKQR ISGNLDSPEG GFDAIMQVAV CGSLIGWRNV TRLLVFSTDA GFHFAGDGKL GGIVLPNDGQ CHLENNMYTM SHYYDYPsia HLVQKLENN IQTIFAVTEE FQPVYKELKN LIPKSAVGTL SANSSNVIQL IIDAYNSLSS EVILENGKLS EGVTSYKSY CKNGVNGTGE NGRKCSNISI GDEVQFEISI TSNKCPKKDS DSKIRPLGF TEEVEVILQY ICECECQSEG IPESPKCHEG NGTFECGACR CNEGRVGRHC ECSTDEVNSE DMDAYCRKEN SSEICSNNGE CVCGQCVCRK RDNTNEIYSG KFCECDNFNCDRSNGLICGG NGVCKCRVCE CNPNYTGSAC DCSLDTSTCE ASNGQICNGR GICECGVCKC TDPKFQGQTC EMCQTCLGVC AEHKECVQCR AFNKGEKKDT CTQECSYFNI TKVESRDKLP QPVQDPVSH CKEKDVDDCW FYFTYSVNGN NEVMVHVVEN PECPTGPDLE HHHHHH
Purity	>95.0% as determined by SDS-PAGE.
Source	Sf9, Baculovirus cells.
Expiration Date	6 months from date of receipt.