

## Product Datasheet

### Human JAM2 (Sf9) protein orb426680

<b>Description</b>	Recombinant of human JAM2 (Sf9) protein
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Target</b>	JAM2 (Sf9)
<b>Form/Appearance</b>	Sterile Filtered colorless solution. JAM2 protein solution (0.5mg/ml) contains 50% glycerol, 1mM DTT & Phosphate Buffered Saline (pH 7.4).
<b>Storage</b>	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
<b>Note</b>	For research use only.
<b>Protein Sequence</b>	ADPFSAPKDQ QVVTAVEYQE AILACKTPKK TVSSRLEWKK LGRSVSFVYY QQTLOGDFKN RAEMIDFNIR IKNVTRSDAG KYRCEVSAPS EQGQNL EEDT VTLEVLVAPA VPSCEVPSSA LSGTVVELRC QDKEGNPAPE YTWFKDGIRL LENPRLGSQS TNSSYTMNTK TGTLOFN TVS KLD TGEYSCE ARNSVGYRRC PGKRMQVDDL NISLEPKSCD KTHTCPPCPA PELLGGPSVF LFPPKPKDTL MISRTPEVTC VVDVSHEDP EVKFNWYVDG VEVHNAKTKP REEQYNSTYR VVSVLTVLHQ DWLNGKEYKC KVSNAKALPAP IEKTISKAKG QPREPQVYTL PPSRDELTKN QVSLTCLVKG FYPSDIAVEW ESNGQPENNY KTTPPVLDS GSFFLYSKLT VDKSRWQQGN VFSCSVMHEA LHNHYTQKSL SLSPGKHHHH HH.
<b>Purity</b>	>90.0% as determined by SDS-PAGE.
<b>Source</b>	Sf9, Baculovirus cells.
<b>Expiration Date</b>	6 months from date of receipt.