

## Product Datasheet

### Human LGALS13 protein orb426856

<b>Description</b>	Recombinant of human LGALS13 protein
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Target</b>	LGALS13
<b>Form/Appearance</b>	Sterile Filtered colorless solution. LGALS13 protein solution (0.5mg/ml) is formulated in 1xPBS buffer pH 7.4.
<b>Storage</b>	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
<b>Note</b>	For research use only.
<b>Protein Sequence</b>	macglvasnl nlkpgceIrv rgevapdaks fvlnlgkdsn nlclhfnprf nahgdantiv cnskdggawg teqreavfpf qpgsvaevci tfdqanltvk lpdgyefkfp nrlnleainy maadgdfkik cvafd.
<b>Purity</b>	>90.0% as determined by SDS-PAGE.
<b>Source</b>	Escherichia Coli.
<b>Expiration Date</b>	6 months from date of receipt.