

Product Datasheet

Malaria Pv. MSP1 protein orb426983

Description	Recombinant of Malaria Pv. MSP1 protein
Conjugation	Unconjugated
Target	Malaria Pv. MSP1
Form/Appearance	Sterile Filtered clear solution. Sterile Filtered solution containing Phosphate buffered saline, 25mM arginine and 1mM EDTA.
Storage	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
Note	For research use only.
Purity	>95% pure as determined by SDS-PAGE.
Expiration Date	6 months from date of receipt.