

Product Datasheet

Human MAPT (383aa) protein orb427015

Description	Recombinant of human MAPT (383aa) protein
Reactivity	Human
Conjugation	Unconjugated
Target	MAPT (383aa)
Form/Appearance	Sterile Filtered colorless solution. The MAPT solution (1 mg/ml) contains 20mM Tris-HCl buffer (pH8.0), 10% glycerol and 0.1M NaCl.
Storage	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
Note	For research use only.
Protein Sequence	MGSSHHHHHH SSGLVPRGSH MAEPRQEFEV MEDHAGTYGL GDRKDQGGYT MHQDQEGDTD AGLKAEEAGI GDTPSLEDEA AGHVTVQARMV SKSKDGTGSD DKKAKGADGK TKIATPRGAA PPGQKGQANA TRIPAKTPPA PKTPPSSGEP PKSGDRSGYS SPGSPGTPGS RSRTPSLPTP PTREPKKVAV VRTPPKSPSS AKSRLQTAPV PMPDLKNVKS KIGSTENLKH QPGGGKVVQII NKKLDLSNVQ SKCGSKDNIK HVPGGGGSVQI VYKPVDSLKV TSKCGSLGNI HHKPGGGQVE VKSEKLDKFD RVQSKIGSLD NITHVPGGGN KKIETHKLT F RENAKAKTDH GAEIVYKSPV VSGDTSRPHL SNVSSTGSID MVDSPQLATL ADEVSASLAK QGL.
Purity	>90.0% as determined by SDS-PAGE.
Source	Escherichia Coli.
Expiration Date	6 months from date of receipt.