

Product Datasheet

MIF protein orb427108

Description	Recombinant of MIF protein
Conjugation	Unconjugated
Target	MIF
Form/Appearance	Sterile Filtered White Lyophilized (freeze-dried) powder. Lyophilized from a 0.2µm filtered concentrated solution in 1
Storage	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
Note	For research use only.
Protein Sequence	MPMFIVNTNV PRASVPEGFL SELTQQLAQA TGKPAQYIAV HVVPDQLMTF SGTSDPCALC SLHSIGKIGG AQNRNYSKLL CGLLSDRLHI SPDRVYINYY DMNAANVGWN GSTFA.
Purity	>95.0% as determined by:(a) Analysis by RP-HPLC.(b) Analysis by SDS-PAGE.
Source	Escherichia Coli.
Solubility (25°C)	It is recommended to reconstitute the lyophilized MIF in sterile 18MΩ-cm H2O not less than 100µg/ml, which can then be further diluted to other aqueous solutions.
Expiration Date	6 months from date of receipt.