

Product Datasheet

Human MYL6 protein orb427237

Description	Recombinant of human MYL6 protein
Reactivity	Human
Conjugation	Unconjugated
Target	MYL6
Form/Appearance	Sterile Filtered colorless solution. The MYL6 solution (0.5 mg/ml) contains 20mM Tris-HCl buffer (pH 8.0), 1mM DTT and 10% glycerol.
Storage	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
Note	For research use only.
Protein Sequence	MGSSHHHHHH SSGLVPRGSH MCDFTEDQTA EFKEAFQLFD RTGDGKILYS QCGDVMRALG QNPTNAEVLK VLG NPKSDEM NVKVLDFEHF LPMLQTVAKN KDQGT YEDYV EGLRVFDKEG NGTVMGAEIR HVLVTLGEEKM TEEVEMLVA GHEDSNGCIN YEELVRMVLN G.
Purity	>90.0% as determined by SDS-PAGE.
Source	Escherichia Coli.
Expiration Date	6 months from date of receipt.