

Product Datasheet

Human NANOG (TAT) protein orb427269

Description	Recombinant of human NANOG (TAT) protein
Reactivity	Human
Conjugation	Unconjugated
Target	NANOG (TAT)
Form/Appearance	Sterile Filtered solution. NANOG protein solution contains PBS and 50mM Arginine.
Storage	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
Note	For research use only.
Protein Sequence	VDPACQSL PCFEASDCKE SSPMPVICGP EENYPSLQMS SAEMPHTETV SPLSSMDLL IQDSPDSSTS PKGKQPTSAE NSVAKKEDKV PVKKQKTRTV FSSTQLCVLN DRFQRQKYL LQQMQELSNL LNLSYKQVKT WFQNQRMKSK RWQKNNWPKN SNGVTQKASA PTYPSLYSSY HQGCLVNPTG NLPMWSNQTW NNSTWSNQTQ NIQSWSNHSW NTQTWCTQSW NNQAWNSPFY NCGEESLQSC MQFQPNPAS DLEAALEAAG EGLNVIQQT RYFSTPQTMD LFLNYSNMNQ PEDVGGYGRK KRRQRR
Purity	>95.0% as determined by SDS-PAGE analysis.
Source	Escherichia Coli.
Expiration Date	6 months from date of receipt.