

## Product Datasheet

### Human NEIL2 protein orb427326

<b>Description</b>	Recombinant of human NEIL2 protein
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Target</b>	NEIL2
<b>Form/Appearance</b>	Sterile filtered colorless solution. The NEIL2 solution (0.5mg/ml) contains 20mM Tris-HCl buffer (pH 8.0), 1mM DTT, 10% glycerol and 0.1M NaCl.
<b>Storage</b>	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
<b>Note</b>	For research use only.
<b>Protein Sequence</b>	MGSSHHHHHH SSGLVPRGSH MGSHPMEGPL VRKFHHLVSP FVGQQVVKTG GSSKKLQPAS LQSLWLQDTQ VHGKKLFLRF DLDEEMGPPG SSPTPEPPQK EVQKEGAADP KQVGEPSPGQK TLDGSSRSAAE LVPQGEDDSE YLERDAPAGD AGRWLRVSFG LFGSVVVNDF SRAKKANKRG DWRDPSPRLV LHFGGGGFLA FYNCQLSWSS SPVVTPTCDI LSEKFHRGQA LEALGQAQPV CYTLLDQRYF SGLGNIKNE ALYRAGIHPL SLGSVLSASR REVLVDHVVE FSTAWLQGKF QGRPQHTQVY QKEQCPAGHQ VMKEAFGPED GLQRLTWWCP QCQPQLSEEP EQCQFS.
<b>Purity</b>	>85.0% as determined by SDS-PAGE.
<b>Source</b>	Escherichia Coli.
<b>Expiration Date</b>	6 months from date of receipt.