

Product Datasheet

Human NGFR protein orb427342

Description	Recombinant of human NGFR protein
Reactivity	Human
Conjugation	Unconjugated
Target	NGFR
Form/Appearance	Sterile Filtered colorless solution. NGFR protein solution (1mg/ml) contains Phosphate Buffered Saline (pH 7.4) and 10% glycerol.
Storage	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
Note	For research use only.
Protein Sequence	KEACPTGLYT HSGECKKACN LGEGVAQPCG ANQTVCEPCL DSVTFSDVVS ATEPCKPCTE CVGLQSMSAP CVEADDAVCR CAYGYYQDET TGRCEACRVC EAGSGLVFSC QDKQNTVCEE CPDGTYSDEA NHVDPCLPCT VCEDTERQLR ECTRWADAEC EEIPGRWITR STPPEGSDST APSTQEPEAP PEQDLIASTV AGVVTVMGS SQPVVTRGTT DNLEHHHHHH.
Purity	>95.0% as determined by SDS-PAGE.
Source	Sf9, Baculovirus cells.
Expiration Date	6 months from date of receipt.