

Product Datasheet

Human NMT1 protein orb427367

Description	Recombinant of human NMT1 protein
Reactivity	Human
Conjugation	Unconjugated
Target	NMT1
Form/Appearance	Sterile Filtered colorless solution. NMT1 protein solution (1mg/ml) containing Phosphate buffered saline (pH7.4).
Storage	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
Note	For research use only.
Protein Sequence	MGSSHHHHHH SSSLVPRGSH MGSMADSEET AVKPPAPPLP QMMEGNGNGH EHCSDCENEE DNSYNRGGLS PANDTGAKKK KKKQKKKKEK GSETDSAQDQ PVKMNSLPAE RIQEIQAIE LFSVGQGPAK TMEEASKRSY QFWDTPQVPK LGEVVNTHGP VEPKDNIRQ EPYTLPQGFT WDALDLGDRG VLKELYLLN ENYVEDDDNM FRFDYSPEFL LWALRPPGWL PQWHCGVRVV SSRKLVGFI AIPANIHIYD TEKKMVEINF LCVHKKLRSK RVAPVLIREI TRRVHLEGIF QAVYTAGVVL PKPVGTCRYW HRSLNPRKLI EVKFSHLRN MTMQRTMKLY RLPETPKTAG LRPMETKDIP VVHQLLTRYL KQFHLTPVMS QEEVEHWFYP QENIIDTFVV ENANGEVTD FLSFYTLPTI MNHPTHKSLK AAYSFYNVHT QTPLLDLMSD ALVLAKMKG FVFNALDLME NKTFLKLF GIGDGNLQYY LYNWKCPMSG AEKVGLVLQ.
Purity	>95.0% as determined by SDS-PAGE.
Source	Escherichia Coli.
Expiration Date	6 months from date of receipt.