

## Product Datasheet

### Human NTS (Sf9) protein orb427426

<b>Description</b>	Recombinant of human NTS (Sf9) protein
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Target</b>	NTS (Sf9)
<b>Form/Appearance</b>	Sterile Filtered clear solution. NTS protein solution (1mg/ml) containing Phosphate Buffered Saline (pH 7.4) and 10% glycerol.
<b>Storage</b>	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
<b>Note</b>	For research use only.
<b>Protein Sequence</b>	ADPSDSEEM KALEADFLTN MHTSKISKAH VPSWKMTLLN VCSLVNNLNS PAEETGEVHE EELVARRKLP TALDGFSLA MLTIYQLHKKI CHSRAFAQHWE LIQEDILDTG NDKNGKEEVI KRKIPYILHH HHHH.
<b>Purity</b>	>90.0% as determined by analysis by SDS-PAGE.
<b>Source</b>	Sf9, Baculovirus cells.
<b>Expiration Date</b>	6 months from date of receipt.