

Product Datasheet

Human OAZ1 protein orb427455

Description	Recombinant of human OAZ1 protein
Reactivity	Human
Conjugation	Unconjugated
Target	OAZ1
Form/Appearance	Sterile Filtered colorless solution. The OAZ1 solution (0.25mg/1ml) contains 20mM Tris-HCl buffer (pH 8.0), 0.15M NaCl, 1mM DTT and 30% glycerol.
Storage	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
Note	For research use only.
Protein Sequence	MGSSHHHHHH SSGLVPRGSH MGSMVKSSLQ RILNSHCFAR EKEGDKPSAT IHASRTMPLL SLHSRGGSSS ESSRVSLHCC SNPGPGPRWC SDAPHPPLKI PGGRGNSQRD HNLSANLFYS DDRLNVTEEL TSNDKTRILN VQSRLTDAKR INWRTVLSSG SLYIEIPGGA LPEGSKDSFA VLLEFAEEQL RADHVFCIFH KNREDRAALL RTFSFLGFEI VRPGHPLVPK RPDACFMAYT FERESSGEEE E
Purity	>90% as determined by SDS-PAGE.
Source	E.coli.
Expiration Date	6 months from date of receipt.