

## Product Datasheet

### Human PDCD5 protein orb427568

<b>Description</b>	Recombinant of human PDCD5 protein
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Target</b>	PDCD5
<b>Form/Appearance</b>	Sterile Filtered colorless solution. PDCD5 protein solution contains 1x PBS pH-7.4.
<b>Storage</b>	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
<b>Note</b>	For research use only.
<b>Protein Sequence</b>	MADEELEALR RQLAELQAK HGDPGDAQQ EAKHRGAEMR NSILAQVLDQ SARARLSNLA LVKPEKTKAV ENYLIQMARY GQLSEKVSEQ GLIEILKKVS QQTEKTTTVK LNRRKVMDS DEDDY.
<b>Purity</b>	>95.0% as determined by SDS-PAGE.
<b>Source</b>	Escherichia Coli.
<b>Expiration Date</b>	6 months from date of receipt.