

## Product Datasheet

### Human PDGFD protein orb427585

<b>Description</b>	Recombinant of human PDGFD protein
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Target</b>	PDGFD
<b>Form/Appearance</b>	Sterile Filtered colorless solution. The PDGFD solution (1mg/1ml) contains 20mM Tris-HCl buffer (pH 8.0), 0.4M Urea and 10% glycerol.
<b>Storage</b>	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
<b>Note</b>	For research use only.
<b>Protein Sequence</b>	MGSSHHHHHH SSGLVPRGSH MGSMSYHDR KSKVDLDRLN DDAKRYSTP RNYSVNIREE LKLANVFFP RCLLVQRGG NCGCGTVNWR SCTCNSGKTV KKYHEVLQFE PGHIKRRGRA KTMALVDIQL DHHERCDCIC SSRPPR.
<b>Purity</b>	>90% as determined by SDS-PAGE.
<b>Source</b>	E.coli.
<b>Expiration Date</b>	6 months from date of receipt.