

Product Datasheet

Human PLAC8 protein orb427727

Description	Recombinant of human PLAC8 protein
Reactivity	Human
Conjugation	Unconjugated
Target	PLAC8
Form/Appearance	Sterile Filtered clear solution. PLAC8 protein solution (1mg/ml) containing 20mM Tris-HCl buffer (pH8.0), 0.4M urea and 10% glycerol.
Storage	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
Note	For research use only.
Protein Sequence	MGSSHHHHHH SSGLVPRGSH MGSMQAQAPV VVVTQPGVGP GPAPQNSNWQ TGMCDGCFSDC GVCLCGTFCF PCLGCQVAAD MNECCLCGTS VAMRTLYRTR YGIPGSICDD YMATLCCPHC TLCQIKRDIN RRRAMRTF.
Purity	>85.0% as determined by SDS-PAGE.
Source	Escherichia Coli.
Expiration Date	6 months from date of receipt.