

Product Datasheet

Human POLR1C protein orb427768

Description	Recombinant of human POLR1C protein
Reactivity	Human
Conjugation	Unconjugated
Target	POLR1C
Form/Appearance	Sterile Filtered colorless solution. POLR1C protein solution (0.5mg/ml) containing 20mM Tris-HCl buffer (pH 8.0), 0.15M NaCl, 20% glycerol and 1mM DTT.
Storage	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
Note	For research use only.
Protein Sequence	MGSSHHHHHH SSGLVPRGSH MGSMAASQAV EEMRSRVVLG EFGVRNVHTT DFPGNYSYD DAWDQDRFEK NFRVDVVHMD ENSLEFDMVG IDAAIANAFR RILLAEVPTM AVEKVLVYNN TSIVQDEILA HRLGLIPIHA DPRLFYRNQ GDEEGTEIDT LQFRLQVRCT RNPAAKDSS DPNELYVNHK VYTRHMTWIP LGNQADLFPE GTIRPVHDDI LIAQLRPGQE IDLLMHCVKG IGKDHAKFSP VATASYRLLP DITLLEPVEG EAAEELSRCF SPGVIEVQEV QGKKVARVAN PRLDTSREI FRNEKLKKVV RLARVRDHYI FSVESTGVLP PDVLVSEAIK VLMGKRRFL DELDAVQMD.
Purity	>80.0% as determined by SDS-PAGE.
Source	Escherichia Coli.
Expiration Date	6 months from date of receipt.