

## Product Datasheet

### Human POLR3F protein orb427779

<b>Description</b>	Recombinant of human POLR3F protein
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Target</b>	POLR3F
<b>Form/Appearance</b>	Sterile Filtered colorless solution. The POLR3F solution (0.5mg/ml) contains 20mM Tris-HCl buffer (pH 8.0), 0.2M NaCl, 30% glycerol and 1mM DTT.
<b>Storage</b>	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
<b>Note</b>	For research use only.
<b>Protein Sequence</b>	MGSSHHHHHH SSGLVPRGSH MGSMAEVKVK VQPPDADPVE IENRIIELCH QFPHGITDQV IQNEMPHIEA QQRAVAINRL LSMGQLDLLR SNTGLLYRIK DSNAGKMKG SDNQEKLVIYQ IIEDAGNKG WSRDIRYKSN LPLTEINKIL KNLESKLIK AVKSVAASKK KVMYLYNLQP DRSVTGGAWY SDQDFESEFV EVLNQOCFKF LQSKAETARE SKQNPMIQRN SSFASSHEVW KYICELGISK VELSMEDIET ILNTLIYDGK VENTIIAAKE GTVGSVDGHM KLYRAVNPII PPTGLVRAPC GLCPVFDDCH EGGEISPSNC IYMTEWLEF.
<b>Purity</b>	>90% as determined by SDS-PAGE.
<b>Source</b>	E.coli.
<b>Expiration Date</b>	6 months from date of receipt.