

Product Datasheet

Human PPIH protein orb427800

Description	Recombinant of human PPIH protein
Reactivity	Human
Conjugation	Unconjugated
Target	PPIH
Form/Appearance	Sterile filtered colorless solution. 1 mg/ml solution containing 1x PBS pH-7.4 10% glycerol.
Storage	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
Note	For research use only.
Protein Sequence	MAVANSSPVN PVVFFDVSIG GOEVGRMKIE LFADVVPKTA ENFRQFCTGEFRKDGVPICY KGSTFHRVIK DFMIQGGDFV NGDGTGVASI YRGPFADEF KLRHSAPGLL SMANSGPSTN GCQFFITCSK CDWLDGKHVV FGKIIDGLLV MRKIENVPTG PNNKPKLPVV ISQCGEM.
Purity	>95.0% as determined by:(a) Analysis by RP-HPLC.(b) Analysis by SDS-PAGE.
Source	Escherichia Coli.
Activity	Specific activity is > 220 nmoles/min/mg, and is defined as the amount of enzyme that cleaves 1umole of suc-AAFP-pNA per minute at 25C in Tris-Hcl pH8.0 using chymotrypsin.
Expiration Date	6 months from date of receipt.