

## Product Datasheet

### Human PSPH protein orb427950

<b>Description</b>	Recombinant of human PSPH protein
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Target</b>	PSPH
<b>Form/Appearance</b>	Sterile filtered colorless solution. The protein contains 20mM Hepes pH 7.5, 1mM DTT &100mM KCl2.
<b>Storage</b>	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
<b>Note</b>	For research use only.
<b>Protein Sequence</b>	MVSHSELRKL FYSADAVCFD VDSTVIREEG IDELAKICGV EDAVSEMTRR AMGGAVPFKA ALTERLALIQ PSREQVQRLLI AEQPPHLTPG IRELVSRMQE RNVQVFLISG GFRSIVEHVA SKLNIPATNV FANRLKFYFN GEYAGFDETQ PTAESGGKGGK VIKLLKEKFH FKKIIMIGDG ATDMEACPPA DAFIGFGGNV IRQVKDNAK WYITDFVELL GELEE.
<b>Purity</b>	>95.0% as determined by:(a) Analysis by RP-HPLC.(b) Analysis by SDS-PAGE.
<b>Source</b>	Escherichia Coli.
<b>Expiration Date</b>	6 months from date of receipt.