


## Product Datasheet

### Human Alpha 1B glycoprotein ELISA Kit orb397030

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Sandwich ELISA for quantitative detection of Human Alpha 1B glycoprotein concentrations.  |
| <b>Reactivity</b>          | Human   |
| <b>Range</b>               | 31.2 pg/ml - 2000 pg/ml   |
| <b>Sensitivity</b>         | 15 pg/ml  |
| <b>Tested Applications</b> | ELISA   |
| <b>Target</b>              | A1BG  |
| <b>Storage</b>             | Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.   |
| <b>Note</b>                | For research use only.  |
| <b>Expiration Date</b>     | Please enquire.   |
| <b>Application Notes</b>   | <b>Detection Method:</b> Colorimetric<br><b>Detection Range:</b> Cell Culture Supernates, Serum, Plasma, Cell Lysates, Tissue Homogenates |



Representative Standard Curve