

Product Datasheet

Human RAB31 protein orb428017

Description	Recombinant of human RAB31 protein
Reactivity	Human
Conjugation	Unconjugated
Target	RAB31
Form/Appearance	Sterile Filtered clear solution. RAB31 protein solution (0.5mg/ml) containing 20mM Tris-HCl buffer (pH 8.0), 0.1M NaCl and 30% glycerol.
Storage	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
Note	For research use only.
Protein Sequence	MRGSHHHHHH GMASMTGGQQ MGRDLYDDDD KDRWGSHMMA IRELKVCLLG DTGVGKSSIV CRFVQDHFDH NISPTIGASF MTKTVP CGNE LHKFLIWDTA GQERFHSLAP MYYRGSAAAV IVYDITKQDS FYTLKKWVKE LKEHGPENIV MAIAGNKCDL SDIREVPLKD AKEYAESIGA IVVETSAKNA INIEELFQGI SRQIPPLDPH ENGNGTIKV EKPTMQASRR CC.
Purity	>90% as determined by SDS-PAGE.
Source	Escherichia Coli.
Expiration Date	6 months from date of receipt.