

## Product Datasheet

### S.Typhi OMP protein orb428246

<b>Description</b>	Recombinant of S.Typhi OMP protein
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA
<b>Target</b>	S.Typhi OMP
<b>Form/Appearance</b>	Sterile Filtered solution. Sterile Filtered solution containing 10mM Tris-HCl, 1mM EDTA and 50mM arginine.
<b>Storage</b>	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
<b>Note</b>	For research use only.
<b>Purity</b>	>95% pure as determined by 12% PAGE (coomassie staining).
<b>Source</b>	Escherichia Coli.
<b>Expiration Date</b>	6 months from date of receipt.