

Product Datasheet

Human SCGB1A1 (His) protein orb428327

Description	Recombinant of human SCGB1A1 (His) protein
Reactivity	Human
Conjugation	Unconjugated
Target	SCGB1A1 (His)
Form/Appearance	Filtered White lyophilized (freeze-dried) powder. Filtered (0.4 µm) and lyophilized from 0.5 mg/ml in 0.05M phosphate buffer and 0.075M NaCl, pH 7.4.
Storage	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
Note	For research use only.
Protein Sequence	MKHHHHHHASEICPSFORVI ETLLMDTPSS YEAAMELFSP DQDMREAGAQ LKKLVDTLPQ KPRESIIKLM EKIAQSSLCN.
Purity	>90.0% as determined by SDS-PAGE.
Source	Escherichia Coli.
Solubility (25°C)	It is recommended to add 200µl deionized water to prepare a working stock solution of approximately 0.5mg/ml and let the lyophilized pellet dissolve completely. SCGB1A1 is not sterile! Please filter the product by an appropriate sterile filter before use.
Expiration Date	6 months from date of receipt.