

Product Datasheet

Human SDSL protein orb428361

Description	Recombinant of human SDSL protein
Reactivity	Human
Conjugation	Unconjugated
Target	SDSL
Form/Appearance	Sterile Filtered colorless solution. The SDSL solution (0.5mg/ml) contains 20mM Tris-HCl buffer (pH 8.0), 2mM DTT, 0.1M NaCl and 10% glycerol.
Storage	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
Note	For research use only.
Protein Sequence	MGSSHHHHHH SSGLVPRGSH MGSMDGPVA EHAKQEPFHV VTPLESWAL SQVAGMPVFL KCENVQPSGS FKIRGIGHFC QEMAKKGRH LVCSSGGNAG IAAAYAARKL GIPATIVLPE STSLQVVQRL QGEGAEVQLT GKVWDEANLR AQELAKRDGW ENVPPFDHPL IWKGHASLVQ ELKAVLRTPP GALVLAVGGG GLLAGVVAGL LEVGWQHVP I AMETHGAHC FNAAITAGKL VTLPDITSVA KSLGAKTVAA RALECMQVCK IHSEVVEDTE AVSAVQQLD DERMLVEPAC GAALAAIYSG LLRRLQAEGC LPPSLTSVVV IVCGGNNINS RELQALKTHL GQV.
Purity	>90% as determined by SDS-PAGE.
Source	E.coli.
Expiration Date	6 months from date of receipt.