


## Product Datasheet

### Human IgG ELISA Kit orb397065

<b>Description</b>	Sandwich ELISA for quantitative detection of Human IgG concentrations.
<b>Reactivity</b>	Human
<b>Range</b>	1560 pg/ml - 100000 pg/ml
<b>Sensitivity</b>	800 pg/ml
<b>Tested Applications</b>	ELISA
<b>Target</b>	IgG
<b>Storage</b>	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
<b>Note</b>	For research use only.
<b>Expiration Date</b>	Please enquire.
<b>Application Notes</b>	<b>Detection Method:</b> Colorimetric <b>Detection Range:</b> Cell Culture Supernates, Serum, Plasma, Cell Lysates, Tissue Homogenates



Representative Standard Curve