

## Product Datasheet

### Human Septin-6 protein orb428378

<b>Description</b>	Recombinant of human Septin-6 protein
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Target</b>	Septin-6
<b>Form/Appearance</b>	Sterile Filtered colorless solution. The SEPT6 solution (0.5mg/ml) contains 20mM Tris-HCl buffer (pH 8.0), 2mM DTT, 0.1M NaCl and 10% glycerol.
<b>Storage</b>	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
<b>Note</b>	For research use only.
<b>Protein Sequence</b>	MGSSHHHHHH SSGLVPRGSH MGSMAATDIA RQVGEGCRTV PLAGHVGFDS LPDQLVNKSV SQGFCFNILC VGETGLGKST LMDTLFNTKF EGEPATHTQP GVQLQSNTYD LQESNVRLKL TIVSTVGFGD QINKEDSYKP IVEFIDAQFE AYLQEELKIR RVLHTYHDSR IHVCLYFIAP TGHSLKSLDL VTMKKLDSKV NIIPAAKAD AISKSELTKE KIKITSELVS NGVQIQFPT DDESVAEING TMNAHLPFVAV IGSTEELKIG NKMMRARQYP WGTVQVENEH HCDFVKLREM LIRVNMEDLR EQHTRHYEL YRRCKLEEMG FKDTDPDSKP FSLQETYEAK RNEFLGELQK KEEEMRQMFV QRVKEKEAEL KEAEKELHEK FDRLKLLHQD EKKKLEDKSK SLDDEVNAFK QRKTAELPQ SQGSQAGGSQ TLKRDKEKKN NPWLCTE
<b>Purity</b>	>90% as determined by SDS-PAGE.
<b>Source</b>	E.coli.
<b>Expiration Date</b>	6 months from date of receipt.