

Product Datasheet

Human SLC4A4 protein orb428475

Description	Recombinant of human SLC4A4 protein
Reactivity	Human
Conjugation	Unconjugated
Target	SLC4A4
Form/Appearance	Sterile filtered colorless solution. The SLC4A4 protein solution contains 20mM Tris-HCl, pH-7.5, 0.5mM DTT, 0.1mM PMSF and 10% glycerol.
Storage	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
Note	For research use only.
Protein Sequence	MRGSHHHHHH GMASMTGGQQ MGRDLYDDDD KDRWGSMS T NVEGKPSNLG ERGRARSSTF LRVVQPMFNH SIFTS AVSPA AERIRFILGEEDDSPAPPQL FTELDELLAV DGQEMEWKET ARWIKFEEKV EQGGERWSKP HVATLSLHSL FELRTCMEKG SIMLDREASS LPQLVEMIVDHQIETGLLKP ELKDKVITYTL LRKHRHQTKK SNLRSLADIG KTVSSASRMF TNPDNGSPAM THRNLTSSSL NDISDKPEKD QLKNKFMKKLPRDAEASNVL VGEVDFLDTP FIAFVRLQQA VMLGALTEVP VPTRFLFILL GPKGKAksyH EIGRAIATLM SDEVFHDIA Y KAKDRHDLIAGIDFLDEVI VLPPGEWDPA IRIEPPKSLP SSDKRKNMYS GGENVQMNGD TPHDGGHGGG GHGDCEELQR TGRFCGGLIK DIKRKAPFFASDFYDALNIQ.
Purity	>90.0% as determined by SDS-PAGE.
Source	Escherichia Coli.
Expiration Date	6 months from date of receipt.