


Product Datasheet

Human PDGF-AA ELISA Kit orb397076

Description	Sandwich ELISA for quantitative detection of Human PDGF-AA concentrations.
Reactivity	Human
Range	15.6 pg/ml - 1000 pg/ml
Sensitivity	7 pg/ml
Tested Applications	ELISA
Target	PDGFA
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	Detection Method: Colorimetric Detection Range: Cell Culture Supernates, Serum, Plasma



Representative Standard Curve