

Product Datasheet

Human SNCA Delta-NAC protein orb428504

Description	Recombinant of human SNCA Delta-NAC protein
Reactivity	Human
Conjugation	Unconjugated
Target	SNCA Delta-NAC
Form/Appearance	Sterile filtered colorless solution. The SNCA Delta-NAC protein solution (1mg/ml) contains 20mM Tris-HCl buffer pH 7.5 and 100mM NaCl.
Storage	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
Note	For research use only.
Protein Sequence	MDVFMKGLSK AKEGVVAAAE KTKQGVAEAA GKTKEGVLYV GSKTKEGVVH GVATVAEKTG GTEIWMKKDQ LGKNEEGAPQ EGILEDMPVD PDNEAYEMPS EEGYQDYEP E A.
Purity	>95.0% as determined by SDS-PAGE.
Source	Escherichia Coli.
Expiration Date	6 months from date of receipt.